



More energy. More adventure.



More family. More time.



Stylish on the Outside.  
Smart on the Inside.

For more information please contact your  
local Fisher & Paykel Healthcare representative

**Manufacturer** Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348, Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com  
**Australia** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22  
**Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **China** Tel: +86 20 3205 3486 Fax: +86 20 3205 2132  
**France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7181 98599 0 Fax: +49 7181 98599 66  
**India** Tel: +91 80 4284 4000 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **Northern Ireland** Tel: +0800 132 189  
**Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10  
**Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625  
**Turkey** Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** (EU Authorised Representative) Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: +44 1628 626 136 Fax: +44 1628 626 146  
**USA** Tel: +1 800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001



www.fphcare.com

**Fisher & Paykel**  
HEALTHCARE

**Fisher & Paykel**  
HEALTHCARE

185045373 REVC © 2010 Fisher & Paykel Healthcare Limited

## Introduce an F&P ICON™ to your family



**ThermoSmart™**  
More humidity.  
More comfort.

**SensAwake™**  
Responsive to  
waking moments.

**InfoSmart™**  
More communication.  
More results.

## Great design enriches your life

Your quality of life shouldn't be compromised by CPAP therapy. Your sleep surroundings needn't be either.

The F&P ICON™ has been designed from the outside-in to answer the CPAP user's strong call for a compact, stylish CPAP that blends seamlessly into any bedroom environment.

Strong on looks and even stronger on intelligence, you can feel at ease knowing you are being taken care of by a full range of the world's leading clinical technologies for treating Sleep Apnea. The F&P ICON™ has a range of useful features that compliment your CPAP therapy as well as everyday life.

### F&P ICON™ Lifestyle Features

<b>Size:</b>	Small and compact for the night-stand, with a fully integrated humidifier, power supply and a hidden water chamber so you still get the best in humidification and comfort while you are sleeping.
<b>Style:</b>	Designed to look more like a household appliance rather than a medical device with neutral colours to blend into any bedroom environment.
<b>Forward Facing Display:</b>	You now have complete control of the menu system at the touch of a button without having to sit up in bed to find the controls.
<b>Clock/Alarm:</b>	No need to squeeze multiple items on the night-stand. The F&P ICON™ is all you need with a clock and alarm function.
<b>Sound:</b>	AlarmTunes™ means that you can wake-up to your favourite songs! Download music onto your SmartStick™ using SmartStick™ Studio software and you will enjoy less 'tone' and more 'tunes'.

## Focus on features

**Water Chamber**  
One-piece chamber. Holds 420 ml of water and is dishwasher safe.

**Chamber Lid**  
Must be correctly fitted for optimal pressure delivery and can also be used as a handle to carry the device.

**Menu System**  
The Home Screen with a clock is the default display. Each point of the clock-face represents the available settings.

**SmartDial™**  
Intuitive controls allow one-touch navigation of the Menu System, using a Press, Turn, Press principle.

**Air Filter**  
Ensures the air is filtered from particles and dust.

**SmartStick™ / USB Port**  
The USB Port houses the SmartStick™ to record patient compliance, efficacy and flow data.

**Elbow**  
Used to connect a ThermoSmart™ or standard Breathing Tube to the device. The Elbow also rotates for a preferred position.